

INTRODUCING WATERBATH SERIES



The Water Bath is designed to take care of heating and cooling applications. It uses microprocessor based technology and is fitted with a powerful pump to take care of both internal and external circulation applications. It incorporates a heater to take care of fast heating rates. The refrigeration system is strong to help achieve fast cooling times. The bath has a drain to help in emptying the circulator tank for periodical cleaning. The grill is removable to help in cleaning the condenser. The Circulators are small and compact and are offered in bath sizes of 3, 5, 9.2 and 20 litres.

FEATURES

- Microprocessor-controlled
- User-friendly keypad
- Warning & safety functions
- Continuous Operation
- External Circulation
- Large LED Display
- Easy to Operate
- Can be used as chiller, circulator, water bath & refrigerated Bath Table Top Model



Specification:

Model	NK5240-V	NK5240-H
Circulator Bath	Vertical 5 Ltr	Horizontal 5 Ltr
Temperature Range	5°C to 85°C	5°C to 85°C
Temperature Stability	+/- 0.3 °C	+/- 0.3 °C
Display	Dual 7 Segment LED	
Display Resolution	0.1 °C	0.1 °C
Temperature Control	PID +	PID +
Heater Capacity	2000 watts	2000 watts
Cooling Capacity @ 19°C	320 watts	320 watts
Pump Flow Rate	12.5 Ltrs/min	12.5 Ltrs/min
Bath Opening [WxLxD cms]	11 x 11 x 15	11 x 11 x 15
Dimensions[WxLxH cms]	42 x 42 x 36	42 x 42 x 15
Refrigerant	R134a	R134a